


UP-D898MD

The Sony UP-D898MD thermal printer is the ideal choice for digital ultrasound systems. This monochrome black and white printer incorporates a USB 2.0 interface for easy integration to digital ultrasound systems. The UP-D898MD printer packs a powerful set of features in a small package. The chassis is extremely small and lightweight with a remarkable print speed of approximately 1.9 seconds.

Digital Graphic Printer



Features:

- Ideal for today's advance ultrasound systems
- User-friendly controls on front panel
- High resolution of 325 DPI
- Compact, space-saving design
- High speed printing of approximately 1.9 seconds
- USB 2.0 Interface
- Utilizes existing 110 Series media by Sony: UPP-110S/5, UPP-110HD, UPP-110HG

B&W**5.8
lb****325
DPI****256
gradation
steps****1.9 secs
per print****Compact
size****USB 2.0
interface**

UP-D898MD Specifications

Dimensions (W x H x D)	6 1/8in x 3 1/2in x 9 1/2in (154mm x 88mm x 240mm)
Weight	5lb 8oz (2.6kg)

Print information

Print technology	Dye Sublimation Thermal Printing
Resolution	325 DPI
Gradations	8 bits (256 levels)
Effective Print Pixels	4096 x 1280 dots
Picture size	12 5/8 x 3 7/8in (Max) 320 x 100mm
Print time	1.9 seconds approximately*
Memory	4096 x 1280 x 8 bit
Printer driver	Windows XP 32/64 bit, Vista 32/64 bit, Windows 7 32/64 bit, Windows 8 32/64 bit
	*A7 print size when using 110HG/110S Media

Interface

USB 2.0 (type A)

Power information

Power requirements	100 to 240V AC, 50/60Hz
Input current	1.3A to 0.6A
(Stand by)	(0.10A to 0.05A)

Operating conditions

Temperature	41°F to 104°F / 5°C to 40°C
Humidity	20% to 80% (noncondensing)

Storage and transport

Temperature	-4°F to +140°F / -20°C to +60°C
Humidity	20% ~ 80%

Standards

Safety	UL 60950-1, 60601-1, ANSI/AAMI ES60601-1, CSA C22.2 No.60950-1, No.601.1, No.60601-1
EMC	FCC Part 15 Subpart-B Digital Device Class A

Accessories

Supplied	(1) Thermal head cleaning sheet (1) CD-ROM (Operating instructions, Printer driver) (1) Service contact list, (1) US Power Cord, (1) USB Cable
----------	--

UP-D898MD Rear Panel



Print Media



UPP-110S/5
Standard Black & White Media

UPP-110HD
High-Density Black & White Media

UPP-110HG
High-Glossy Black & White Media



SYNERGY medical inc

8069 mexico rd
st. peters, mo 63376

office: 636.928.2800

email: info@synergymedco.com

www.synergymedco.com

©2014 Sony Electronics Inc. All rights reserved. Reproduction in whole or in part without written permission is prohibited. All non-metric weights and measurements are approximate. Features and specifications are subject to change without notice. Sony is a trademarks of Sony. All other trademarks are trademarks of their respective owners.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or other appropriately licensed medical professional.
Caution: See product labeling for indications, contraindications, warnings, cautions and directions for use.

ME-D898MD