

MICROFLEX®

NEOGARD®




GREEN
NEOPRENE
NON-STERILE
AMBIDEXTROUS
TEXTURED FINGERS
BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX

POWDER-FREE EXAMINATION GLOVES

Our new Neogard® glove is made from a unique, non latex formulation that delivers outstanding grip in wet and dry environments. With unparalleled tactile sensitivity and comfort, Neogard® gloves are ideal for every task.

- Soft and stretchy chloroprene formulation is highly comfortable with exceptional tactile sensitivity
- Advanced formulation delivers excellent wet and dry grip
- Polymer coated for easy donning, especially on damp hands



A FIRM GRASP ON WET GRIP TECHNOLOGY

PRODUCT CODE		PACKAGING	
XS	C520	100	Gloves/Box
S	C521	10	Boxes/Case
M	C522		
L	C523		
XL	C524		

GLOVE PROFILE		ASTM D 6977	NEOGARD®
Length	mm (inches)	220 (8.7)	245 (9.6)
Cuff Thickness	mm (mils)	n/a	0.07 (2.8)
Palm Thickness	mm (mils)	min. 2.0	0.10 (3.9)
Finger Thickness	mm (mils)	min. 2.0	0.13 (5.1)
Tensile Strength	Before Aging	min. 14	18
	After Aging	min. 14	14
Elasticity/Elongation	Before Aging	min. 500	650
	After Aging	min. 400	500
Pinhole Defect Rate*	(AQL)	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED INDUSTRIES