

MICROFLEX®

▲ COLORTOUCH® PINK



PINK
LATEX
NON-STERILE
AMBIDEXTROUS
FULLY TEXTURED
BEADED CUFF
FOR SINGLE USE ONLY

POWDER-FREE EXAMINATION GLOVES

ColorTouch® Pink gloves provide more than great hand protection; \$1 from every case sold is donated to The Breast Cancer Research Foundation®. These high quality, polymer-coated latex gloves provide outstanding protection and are textured to provide a more secure grip. ColorTouch® Pink gloves slide on easily and are designed to provide exceptional comfort to keep hands relaxed.



- Unique pink color shows your support for a great cause
- Polymer-coated for easy donning
- Fully textured and relaxed comfort for reduced hand fatigue

| PRODUCT CODE | | PACKAGING | |
|--------------|------------|-----------|------------|
| XS | CTP-233-XS | 100 | Gloves/Box |
| S | CTP-233-S | 10 | Boxes/Case |
| M | CTP-233-M | | |
| L | CTP-233-L | | |

| GLOVE PROFILE | | ASTM D 3578 | COLORTOUCH® PINK |
|-----------------------|--------------|-----------------|------------------|
| Length | mm (inches) | min. 220 (8.7) | 245 (9.6) |
| Cuff Thickness | mm (mils) | n/a | 0.10 (3.9) |
| Palm Thickness | mm (mils) | min. 0.08 (3.1) | 0.15 (5.9) |
| Finger Thickness | mm (mils) | min. 0.08 (3.1) | 0.17 (6.7) |
| Tensile Strength | Before Aging | min. 18 | 22 |
| | After Aging | min. 14 | 21 |
| Elasticity/Elongation | Before Aging | min. 650 | 750 |
| | After Aging | min. 500 | 600 |
| Pinhole Defect Rate* | (AQL) | max. 2.5 | 1.5 |

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED INDUSTRIES

Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.

Microflex is donating \$1.00 per case sold to The Breast Cancer Research Foundation®. BCRF is dedicated to preventing breast cancer and finding a cure in our lifetime by funding clinical and translational research worldwide. For more information about BCRF, visit www.bcrfcure.org

® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2014 All Rights Reserved.

COLL-SPEC-CTP233_v1

Ansell