

MICROFLEX®

▲ FREEFORM® EC



BLUE
NITRILE
NON-STERILE
AMBIDEXTROUS
TEXTURED FINGERS
EXTENDED BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX

POWDER-FREE EXTENDED CUFF EXAMINATION GLOVES

These extra strong gloves meet or exceed NFPA standards for emergency medical equipment. The long cuff provides durable, extended coverage for high-risk conditions. Ideal if you need your hands to be comfortable and free to move as you work.

- Durable, extended cuff for coverage in high-risk conditions
- Textured fingertips for enhanced grip



NON-STICK PROPERTIES RESIST TAPE AND ADHESIVES

PRODUCT CODE		PACKAGING	
S	FFE-775-S	50	Gloves/Box
M	FFE-775-M	10	Boxes/Case
L	FFE-775-L		
XL	FFE-775-XL		
XXL	FFE-775-XXL		

GLOVE PROFILE		ASTM D 6319	FREEFORM® EC
Length	<i>mm (inches)</i>	min. 220 (8.7)	290 (11.4)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.09 (3.5)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.12 (4.7)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.16 (6.3)
Tensile Strength	<i>Before Aging</i>	min. 14	31
	<i>After Aging</i>	min. 14	27
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	500
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED INDUSTRIES