

# MICROFLEX®

## ▲ MIDKNIGHT®



BLACK  
NITRILE  
NON-STERILE  
AMBIDEXTROUS  
FULLY TEXTURED  
BEADED CUFF  
FOR SINGLE USE ONLY  
NON-LATEX

### POWDER-FREE EXAMINATION GLOVES

MidKnight® gloves are designed with a distinct color that masks stains and provides improved framing effect. This is the glove of choice for individuals who demand both quality and protection against day-to-day hazards.

- Distinctive black color to mask stains
- Excellent choice for wet environments
- Fully textured for firm and consistent grip



### A LEADER IN LABORATORY AND AUTOMOTIVE MARKETS

PRODUCT CODE		PACKAGING	
XS	MK-296-XS	100	Gloves/Box
S	MK-296-S	10	Boxes/Case
M	MK-296-M		
L	MK-296-L		
XL	MK-296-XL		
XXL	MK-296-XXL		

GLOVE PROFILE		ASTM D 6319	MIDKNIGHT®
Length	mm (inches)	min. 220 (8.7)	245 (9.6)
Cuff Thickness	mm (mils)	n/a	0.08 (3.1)
Palm Thickness	mm (mils)	min. 0.05 (2.0)	0.12 (4.7)
Finger Thickness	mm (mils)	min. 0.05 (2.0)	0.14 (5.5)
Tensile Strength	Before Aging	min. 14	27
	After Aging	min. 14	24
Elasticity/Elongation	Before Aging	min. 500	500
	After Aging	min. 400	400
Pinhole Defect Rate*	(AQL)	max. 2.5	1.5

\* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED  
INDUSTRIES