

MICROFLEX®

▲ SUPRENO® SE



BLUE
 NITRILE
 NON-STERILE
 AMBIDEXTROUS
 TEXTURED FINGERS
 BEADED CUFF
 FOR SINGLE USE ONLY
 NON-LATEX

POWDER-FREE EXTENDED CUFF EXAMINATION GLOVES

Supreno® SE gloves provide the durability needed for demanding applications. Offering all the benefits of a nitrile glove, Supreno® SE meets your requirement of strength and comfort.

- Lowest pinhole rates (AQL 0.65) sets a new standard for protection
- Provides durable, heavy-duty coverage
- Non-stick properties resist tape and adhesives



PRODUCT CODE		PACKAGING	
XS	SU-690-XS	100	Gloves/Box
S	SU-690-S	10	Boxes/Case
M	SU-690-M		
L	SU-690-L		
XL	SU-690-XL		

GLOVE PROFILE		ASTM D 6319	SUPRENO® SE
Length	<i>mm (inches)</i>	min. 220 (8.7)	245 (9.6)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.09 (3.5)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.11 (4.3)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.18 (7.1)
Tensile Strength	<i>Before Aging</i>	min. 14	30
	<i>After Aging</i>	min. 14	30
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	400
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	0.65

* Pinhole Defect Rate is FDA minimum requirement.

