

MICROFLEX®

TRANQUILITY®



WHITE
NITRILE
NON-STERILE
AMBIDEXTROUS
TEXTURED FINGERS
BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX

POWDER-FREE EXAMINATION GLOVES

Tranquility® gloves are designed to deliver the perfect balance between protection and comfort. Along with superior strength, these nitrile gloves conform to the hand allowing for excellent tactile sensitivity and a secure grip. A compact storage box reduces space requirements in the office or lab.



- Compact box for reduced storage space
- Textured fingertips
- Non-stick properties resist tape and adhesives



A LEADER IN DENTAL AND HEALTHCARE MARKETS

PRODUCT CODE		PACKAGING	
XS	TQ-601-XS	100	Gloves/Box
S	TQ-601-S	10	Boxes/Case
M	TQ-601-M		
L	TQ-601-L		
XL	TQ-601-XL		

GLOVE PROFILE		ASTM D 6319	TRANQUILITY®
Length	<i>mm (inches)</i>	min. 220 (8.7)	240 (9.5)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.055 (2.2)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.065 (2.6)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.095 (3.7)
Tensile Strength	<i>Before Aging</i>	min. 14	24
	<i>After Aging</i>	min. 14	26
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	400
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED
INDUSTRIES