

MICROFLEX®

▲ ULTRAFORM®



BLUE
NITRILE
NON-STERILE
AMBIDEXTROUS
TEXTURED FINGERTIPS
BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX

POWDER-FREE EXAMINATION GLOVES

Our unique, soft formulation forms to your hands effortlessly, as if it is a second skin. In field tests, users wearing Ultraform® experienced an increase in fine motor task performance compared to a bare hand.

- Packaged 300 gloves per box to save space and reduce waste
- Textured fingertips for an improved grip



CERTIFIED ERGONOMIC EXAM GLOVE

PRODUCT CODE		PACKAGING		
XS	UF-524-XS	300 ⁺	Gloves/Box	+ XL contains 250 gloves/box
S	UF-524-S	10	Boxes/Case	
M	UF-524-M			
L	UF-524-L			
XL	UF-524-XL			

GLOVE PROFILE		ASTM D 6319	ULTRAFORM®
Length	<i>mm (inches)</i>	min. 220 (8.7)	230 (9.1)
Cuff Thickness	<i>mm (mils)</i>	n/a	0.055 (2.2)
Palm Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.06 (2.4)
Finger Thickness	<i>mm (mils)</i>	min. 0.05 (2.0)	0.08 (3.1)
Tensile Strength	<i>Before Aging</i>	min. 14	25
	<i>After Aging</i>	min. 14	27
Elasticity/Elongation	<i>Before Aging</i>	min. 500	500
	<i>After Aging</i>	min. 400	400
Pinhole Defect Rate*	<i>(AQL)</i>	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED INDUSTRIES