

MICROFLEX®

▲ E-GRIP®



NATURAL
LATEX
NON-STERILE
AMBIDEXTROUS
TEXTURED FINGERS
BEADED CUFF
FOR SINGLE USE ONLY

POWDER-FREE EXAMINATION GLOVES

The micro-textured satin finish provides enhanced gripping power and excellent tactile sensitivity. Beaded cuff for easy donning. 510K approved for medical use.

- Offers enhanced tactile sensitivity and grip



ONE OF OUR MOST POPULAR LATEX GLOVES

PRODUCT CODE		PACKAGING	
XS	L970	100	Gloves/Box
S	L971	10	Boxes/Case
M	L972		
L	L973		
XL	L974		

GLOVE PROFILE		ASTM D 3578	E-GRIP®
Length	mm (inches)	min. 220 (8.7)	245 (9.6)
Cuff Thickness	mm (mils)	n/a	0.10 (3.9)
Palm Thickness	mm (mils)	min. 0.08 (3.1)	0.13 (5.1)
Finger Thickness	mm (mils)	min. 0.08 (3.1)	0.14 (5.5)
Tensile Strength	Before Aging	min. 18	22
	After Aging	min. 14	14
Elasticity/Elongation	Before Aging	min. 650	700
	After Aging	min. 500	600
Pinhole Defect Rate*	(AQL)	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED
INDUSTRIES

Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.