

MICROFLEX®

▲ NITRILE EXAM GLOVES



POWDER-FREE EXAMINATION GLOVES

The soft and stretchy feel of these premium gloves gives superior fit and performance. The latex-free construction offers protection without the risk of allergic reactions associated with latex proteins. Fully textured finish offers an excellent grip.

- Superior puncture and chemical resistance



BLUE
NITRILE
NON-STERILE
AMBIDEXTROUS
TEXTURED FINGERS
BEADED CUFF
FOR SINGLE USE ONLY
NON-LATEX



PRODUCT CODE		PACKAGING	
XS	N850	100	Gloves/Box
S	N851	10	Boxes/Case
M	N852		
L	N853		
XL	N854		
XXL	N855		

GLOVE PROFILE		ASTM D 6319	N85
Length	mm (inches)	min. 220 (8.7)	240 (9.5)
Cuff Thickness	mm (mils)	n/a	0.06 (2.4)
Palm Thickness	mm (mils)	min. 0.05 (2.0)	0.08 (3.1)
Finger Thickness	mm (mils)	min. 0.05 (2.0)	0.11 (4.3)
Tensile Strength	Before Aging	min. 14	min. 16
	After Aging	min. 14	min. 14
Elasticity/Elongation	Before Aging	min. 500	500
	After Aging	min. 400	400
Pinhole Defect Rate*	(AQL)	max. 2.5	1.5

* Pinhole Defect Rate is FDA minimum requirement.



RECOMMENDED
INDUSTRIES