



31.5" 8MP Diagnostic Monitor **32HL512D** (future release)

Monitor is DICOM out-of-the box calibrated with brightness stabilization to ensure confident reading over time. Operates in single screen 8MP or in Picture by Picture mode to meet PACS application display requirements.

- Native Resolution: 3840 x 2160 (16:9)
- Picture by Picture: 2X 1920x2160
- Brightness 450 cd/m² (typical)
- DICOM calibrated: 350 cd/m²
- Brightness Stabilization
- Inputs: DP (x2), HDMI (x1), USB3.0 (x2)
- IEC60601-1 / IEC60601-1-2
- FDA Class II (pending)